

STANDARD DIMENSIONS

PART NUMBER	NOMINAL VOLUME (mL)	OVERFLOW (BRIM) CAPACITY (mL)	A IN (MM)	B IN (MM)	ØC IN (MM)	
150-01-0100	100	165	4.27" (108.5)	2.12" (53.8)	2.19" (55.6)	FIGURE 1
150-01-0250	250	380	5.75" (146.1)	2.12" (53.8)	2.75" (69.9)	
150-01-0500	500	584	6.33" (160.9)	2.12" (53.8)	3.35" (85.1)	
150-01-1000	1000	1200	8.60" (218.4)	2.12" (53.8)	3.98" (101.1)	
150-01-2000	2000	2222	9.77" (248.1)	2.12" (53.8)	5.01" (127.3)	
150-02-0100	100	162	4.27" (108.5)	2.12" (53.8)	2.19" (55.6)	
150-02-0250	250	352	5.75" (146.1)	2.12" (53.8)	2.75" (69.9)	
150-02-0500	500	588	6.33" (160.9)	2.12" (53.8)	3.35" (85.1)	
150-02-1000	1000	1182	8.60" (218.4)	2.12" (53.8)	3.98" (101.1)	
150-02-2000	2000	2240	9.73" (247.0)	2.12" (53.8)	5.05" (127.3)	
150-01-0050	50	62	3.30" (83.8)	1.44" (36.6)	1.44" (36.6)	FIGURE 2
150-02-0050	50	62	3.30" (83.8)	1.44" (36.6)	1.44" (36.6)	

