

STANDARD DIMENSIONS

PART NUMBER	VOLUME	A INCHES (MM)	B INCHES (MM)	C INCHES (MM)	D INCHES (MM)	E INCHES (MM)	F INCHES (MM)	
210-005-21	5 mL	-	0.47" (11.9)	-	2.95" (74.9)	0.39" (9.9)	2.90" (73.7)	FIGURE 1
210-014-21	14 mL	-	0.63" (16.0)	-	3.93" (99.8)	0.54" (13.7)	3.87" (98.2)	
210-015-20	15 mL	0.84" (21.3)	0.63" (16.0)	0.51" (13.0)	4.45" (113.0)	0.54" (13.7)	4.28" (108.7)	FIGURE 2
210-015-30	15 mL	0.84" (21.3)	0.63" (16.0)	0.51" (13.0)	5.09" (129.3)	0.54" (13.7)	4.90" (124.5)	
210-026-20	26 mL	1.08" (27.4)	0.85" (21.6)	0.54" (13.7)	4.57" (116.1)	0.67" (16.9)	4.39" (111.4)	
210-032-20	32 mL	1.03" (26.2)	0.79" (20.1)	0.52" (13.2)	6.04" (153.4)	0.63" (16.1)	5.81" (147.5)	
210-050-20	50 mL	1.42" (36.1)	1.12" (28.4)	0.51" (13.0)	4.83" (122.7)	1.06" (27.0)	4.63" (117.7)	
210-050-30	50 mL	1.42" (36.1)	1.12" (28.4)	0.51" (13.0)	4.70" (119.4)	1.06" (27.0)	4.46" (113.2)	
210-060-20	60 mL	1.24" (31.5)	1.02" (25.9)	0.60" (15.2)	6.64" (168.7)	0.68" (17.2)	5.82" (147.9)	

***NOTE: TUBE IS SHOWN WITH SAVILLEX CLOSURE FOR REFERENCE**

FIGURE 1




FIGURE 2

