

STANDARD DIMENSIONS


PART NUMBER	VOLUME	A INCHES (MM)	B INCHES (MM)	C INCHES (MM)	D INCHES (MM)	E INCHES (MM)	F INCHES (MM)	
200-902-40	0.2 mL	0.49" (12.4)	0.45" (11.4)	0.41" (10.4)	1.36" (34.5)	13.6" (34.5)	2.17" (5.5)	FIGURE 1
200-905-42	0.5 mL	0.49" (12.4)	0.45" (11.4)	0.41" (10.4)	1.35" (34.3)	13.5" (34.3)	2.17" (5.5)	
200-913-44	1.3 mL	0.53" (13.5)	0.48" (12.2)	0.60" (15.2)	1.55" (39.4)	15.5" (39.4)	3.35" (8.5)	
200-004-13	4 mL	0.73" (18.5)	0.59" (15.0)	0.24" (6.1)	2.06" (52.3)	20.6" (52.3)	4.69" (11.9)	
200-013-14	13 mL	1.08" (27.4)	0.87" (22.1)	0.54" (13.7)	2.21" (56.1)	22.1" (56.1)	7.44" (18.9)	FIGURE 2
200-915-44	1.5 mL	0.64" (16.3)	0.50" (12.7)	0.43" (10.9)	1.24" (31.5)	12.4" (31.5)	4.33" (11.0)	
200-002-44	2 mL	0.64" (16.3)	0.48" (12.2)	0.43" (10.9)	1.72" (43.7)	17.2" (43.7)	4.49" (11.4)	
200-004-44	4 mL	0.64" (16.3)	0.56" (14.2)	0.43" (10.9)	2.37" (60.2)	23.7" (60.2)	4.65" (11.8)	FIGURE 3
200-902-41	0.2 mL	0.43" (10.8)	0.45" (11.5)	0.15" (3.8)	1.28" (32.4)	12.8" (32.4)	2.17" (5.5)	
200-905-43	0.5 mL	0.43" (10.8)	0.45" (11.5)	0.15" (3.8)	1.28" (32.4)	12.8" (32.4)	2.17" (5.5)	

***NOTE: AUTOSAMPLER VIALS ARE SHOWN WITH SAVILLEX CLOSURES FOR REFERENCE**

FIGURE 1

FIGURE 2

FIGURE 3



STANDARD DIMENSIONS

PART NUMBER	VOLUME	A INCHES (MM)	B INCHES (MM)	C INCHES (MM)	D INCHES (MM)
200-010-13	10 mL	0.60" (15.2)	0.58" (14.7)	0.09" (2.3)	3.99" (101.3)
200-020-13	20 mL	1.07" (27.2)	1.05" (26.7)	0.08" (2.0)	2.35" (59.7)
200-030-13	30 mL	1.06" (26.9)	1.05" (26.7)	0.08" (2.0)	3.99" (101.3)

