

STANDARD DIMENSIONS

PART NUMBER	VOLUME	A INCHES (MM)	B INCHES (MM)	C INCHES (MM)	D INCHES (MM)	E INCHES (MM)	F INCHES (MM)
210-004-20	4 mL	0.96" (24.4)	0.60" (15.2)	0.40" (10.2)	1.82" (46.2)	0.56" (14.2)	1.63" (41.3)
210-010-30	10 mL	1.44" (36.6)	1.03" (26.2)	0.62" (15.7)	1.60" (40.6)	0.99" (25.1)	1.34" (34.1)
210-017-20	17 mL	1.44" (36.6)	1.03" (26.2)	0.62" (15.7)	1.60" (40.6)	1.01" (25.6)	1.31" (33.4)
210-025-30	25 mL	1.44" (36.6)	1.03" (26.2)	0.06" (1.5)	3.09" (78.5)	0.99" (25.1)	2.74" (69.7)
210-055-11	55 mL	1.31" (33.3)	1.19" (30.2)	0.62" (15.7)	4.31" (109.5)	-	-
210-060-70	60 mL	1.44" (36.6)	1.16" (29.5)	0.62" (15.7)	5.22" (132.6)	-	-
210-075-11	75 mL	-	1.36" (34.5)	-	4.11" (104.4)	1.30" (33.0)	3.71" (94.3)

***NOTE: BLOCK DIGESTION TUBE IS SHOWN WITH SAVILLEX CLOSURE FOR REFERENCE**

