

### STANDARD DIMENSIONS

PART NUMBER	VOLUME	A INCHES (MM)	B INCHES (MM)	C INCHES (MM)	D INCHES (MM)	E INCHES (MM)	F INCHES (MM)
200-702-42	0.2 mL	0.49" (12.4)	0.42" (10.7)	0.41" (10.4)	1.36" (34.5)	2.05" (5.2)	6.93" (17.6)
200-705-42	0.5 mL	0.49" (12.4)	0.42" (10.7)	0.41" (10.4)	1.36" (34.5)	2.05" (5.2)	12.0" (30.6)

**\*NOTE: FORENSIC VIALS ARE SHOWN WITH SAVILLEX CLOSURES FOR REFERENCE**

