


STANDARD DIMENSIONS

PART NUMBER	VOLUME	A INCHES (MM)	B INCHES (MM)	C INCHES (MM)	D INCHES (MM)	E INCHES (MM)	F INCHES (MM)
201-060-12-033-26	60 mL	1.44" (36.6)	1.21" (30.7)	6.85" (174.0)	1.31" (33.3)	4.92" (125.0)	1.00" (25.3)
201-060-12-033-27	60 mL	1.47" (37.3)	1.21" (30.7)	6.99" (174.6)	1.46" (37.1)	4.92" (125.0)	1.00" (25.3)
201-060-12-033-36	60 mL	1.44" (36.6)	1.21" (30.7)	6.85" (174.0)	0.57" (14.5)	4.92" (125.0)	1.00" (25.3)
201-060-12-033-37	60 mL	1.47" (37.3)	1.21" (30.7)	6.83" (173.5)	0.57" (14.5)	4.92" (125.0)	1.00" (25.3)
501-275-02-058-42	375 mL	2.29" (58.2)	1.94" (49.3)	12.99" (330.0)	1.10" (28.0)	11.49" (291.8)	1.60" (40.6)

***NOTE: IMPINGERS ARE SHOWN WITH SAVILLEX CLOSURES FOR REFERENCE**



TITLE: **IMPINGER DIMENSIONS**