

STANDARD DIMENSIONS

PART NUMBER	VOLUME	A INCHES (MM)	B INCHES (MM)	C INCHES (MM)	D INCHES (MM)	E INCHES (MM)	F INCHES (MM)
200-905-50	0.5 mL	0.40" (10.0)	0.31" (7.8)	0.08" (2.0)	1.22" (30.9)	2.44" (6.2)	10.3" (26.2)
200-915-50	1.5 mL	0.64" (16.3)	0.50" (12.7)	0.43" (10.9)	1.67" (42.3)	3.35" (8.5)	14.8" (37.5)

***NOTE: MICRO VIAL IS SHOWN WITH SAVILLEX CLOSURE FOR REFERENCE**

