

STANDARD DIMENSIONS

	PART NUMBER	*PSI RATING 100°F/36°C	VOLUME	A INCHES (MM)	B INCHES (MM)
FIGURE 1	102-0010-01	100 PSI	10 ml	3.06" (77.7)	1.02" (25.9)
	100-0125-01	50 PSI	125 ml	5.06" (128.5)	1.97" (50.0)
FIGURE 2	102-0475-01	25 PSI	475 ml	4.55" (115.6)	3.35" (85.1)
	102-0475-03	25 PSI	475 ml	4.55" (115.6)	3.35" (85.1)
	102-1100-01	75 PSI	1100 ml	6.62" (168.1)	4.22" (107.2)
	102-1100-03	75 PSI	1100 ml	6.45" (163.8)	4.24" (107.7)

***NOTE: FOR EVERY 2°F/1°C INCREASE ABOVE 100°F/38°C, REDUCE PRESSURE RATING BY 1 PSI**

FIGURE 1





FIGURE 2



SavilleX 6133 Baker Road
 Minnetonka, MN
 55345-5910 USA
 (952)935-4100

TITLE:
**SEALED
 JAR DIMENSIONS**