

STANDARD DIMENSIONS

PART NUMBER	VOLUME	A INCHES (MM)	B INCHES (MM)	C INCHES (MM)	D INCHES (MM)
210-004-20	4 ml	0.96" (24.4)	0.60" (15.2)	0.40" (10.2)	1.82" (46.2)
210-010-30	10 ml	1.44" (36.6)	1.03" (26.2)	0.62" (15.7)	1.60" (40.6)
210-017-20	17 ml	1.44" (36.6)	1.03" (26.2)	0.62" (15.7)	1.60" (40.6)
210-025-30	25 ml	1.44" (36.6)	1.03" (26.2)	0.62" (15.7)	3.09" (78.5)
210-055-11	55 ml	1.31" (33.3)	1.19" (30.2)	0.06" (1.5)	4.31" (109.5)
210-060-70	60 ml	1.44" (36.6)	1.16" (29.5)	0.62" (15.7)	5.22" (132.6)
210-075-11	75 ml	1.44" (36.6)	1.36" (34.5)	0.23" (5.8)	4.11" (104.4)

***NOTE: BLOCK DIGESTION TUBE IS SHOWN WITH SAVILLEX CLOSURE FOR REFERENCE**

