

### STANDARD DIMENSIONS

PART NUMBER	PSI RATING 73°F/23°C	A INCHES (MM)	B INCHES (MM)
600-023-01	60 PSI	0.40" (10.2)	0.96" (24.3)
600-024-01	25 PSI	0.54" (13.7)	1.08" (27.4)
600-033-01	10 PSI	0.62" (15.7)	1.44" (36.6)
600-033-71	10 PSI	0.62" (15.7)	1.44" (36.6)

**NOTES:**

1. CONSULT SAVILLEX FOR ELEVATED TEMPERATURE RATINGS.
2. PSI RATINGS ARE BASED ON ABOVE CLOSURES WITH MATING SAVILLEX VESSELS

