


STANDARD DIMENSIONS

PART NUMBER	A INCHES (MM)	B INCHES (MM)	C INCHES (MM)	D INCHES (MM)
730-2001	0.70" (17.8)	6.88" (174.8)	1.00" (25.4)	2.71" (68.7)
730-2002	1.27" (32.3)	8.57" (217.7)	1.39" (35.3)	2.71" (68.7)




 10321 West 70th St.
 Eden Prairie, MN
 55344-3446 USA
 (952)935-4100

TITLE: **VIAL TRAY DIMENSIONS**